

**Product Name :**  
Hydraulic Jacks

**Product Code :**  
NLAB-MACHINELABWK28004



**Description :**

Hydraulic Jacks

**Technical Specification :**

Hydraulic Jacks .

~~250~~ Capacity  
in Tons

~~250~~ Maximum  
(~~200~~)

~~150~~ Hydraulic  
Lift

~~300~~ew  
Extn.  
(-)

~~600~~ Ht.  
(~~600~~)

~~600~~ Dia.

~~600~~ly Dia.

**Size**  
380  
x 126  
**Weight in**  
**Kgs.**  
  
(approx.)

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003