

**Product Name :**  
 Hydraulic Pipe Bending Machine

**Product Code :**  
 NLAB-VOCATIONALTVET13027



**Description :**

Hydraulic Pipe Bending Machine

**Technical Specification :**

Capacity

|                          |  |
|--------------------------|--|
| 29 Ton.                  |  |
| 250 mm Stroke            |  |
| 210, 225 mm (mm.)        |  |
| Bending Formers (D, D)½, |  |

|                      |  |
|----------------------|--|
| 3, 4                 |  |
| 275 mm x 6.5 mm Pipe |  |
| 530 x 203            |  |
| 100 kg (Approx.)     |  |

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003