

Product Name :
 Hydraulic Pipe Bending Machine

Product Code :
 NLAB-VOCATIONALTVET13027



Description :

Hydraulic Pipe Bending Machine

Technical Specification :

Capacity

29 Ton.	
250 mm Stroke	
210, 225 mm (mm.)	
Bending Formers (D, D)½,	

3, 4	
275 mm x 6.5 mm Pipe	
530 x 203	
100 kgs (Approx.)	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003