

Product Name :
Hydraulic Pipe Bending Machine

Product Code :
NLAB-MACHINELABWK34017



Description :

Hydraulic Pipe Bending Machine

Technical Specification :

Type: Hydraulic

Max. pressure: 10 ton

Max. oil tank pressure: 42mpa

Max. outside pipe diameter: 2"

Min. outside pipe diameter: ½"

Max. ram stroke: 155mm

Standard Welding Equipment Laboratory Set Hydraulic Pipe Bending Machine
Manufacturer,Supplier and Hydraulic Pipe Bending Machine Exporter in
India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003