

Product Name :
Hydraulic Puller

Product Code :
NLAB-VOCATIONALTVET20012



Description :

Hydraulic Puller

Technical Specification :

Hydraulic Puller Stroke-50 mm, Spread-200 mm, Reach-140 mm.

Compact - The Self Contained Hydraulic Puller - Pump and Cylinder Save Space, As No Separate Pump or Hose is Required.

360° Rotating Pump Provides Flexibility For Using the Puller in the Most Convenient Position.

Piston Rod - Heat Treated and Hard Chrome Plated.

Effortless and Simple Release Valve.

Heat Treated Forged Alloy Steel Jaws.

Fast and Smooth Spring Retraction.

Safety - The Integral Safety Release Valve Prevents Overloading.

Practical - Easy to Convert to 2 Arm or 3 Arm Puller (For Uniform Pulling)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003