

Product Name :
Hydraulic Pullers

Product Code :
NLAB-MACHINELABWK28003



Description :

Hydraulic Pullers

Technical Specification :

Hydraulic Pullers .

Hydraulic puller - Remote controlled comprising of hollow ram cylinder hand operated high pressure, single acting, single piston, Hydraulic pump with pressure gauge adoptor, 2 meter long high pressure rubber hose but without any attachmnets.

Shaft dia. mm

~~20~~ 20

~~25~~ 25-50

~~30~~ 30-90

~~50~~ 50-140

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003