

Product Name :
Hydraulic Shrinkage Mould 25x25x285 mm

Product Code :
NLAB-CEMETLAB130016



Description :

Hydraulic Shrinkage Mould 25x25x285 mm

Technical Specification :

This mould is the same used for the Expansion of Portland cement with the High pressure autoclave, conforming to ASTM C151 and ASTM C490.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003