

Product Name :
Hydraulic Shrinkage Mould 40x40x160 mm

Product Code :
NLAB-CEMETLAB130015



Description :

Hydraulic Shrinkage Mould 40x40x160 mm

Technical Specification :

Three gang mould, individually tested with certified instruments. Made of special alloy steel, minimum hardness HV200. The shrinkage measurement has to be done with the Length comparator with the .This version, conforming to UNI 6687, slightly differ from other version and it is used for the same determination.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003