

Product Name :
Hydraulic Shrinkage Mould 40x40x160 mm

Product Code :
NLAB-CEMETLAB130015



Description :

Hydraulic Shrinkage Mould 40x40x160 mm

Technical Specification :

Three gang mould, individually tested with certified instruments. Made of special alloy steel, minimum hardness HV200. The shrinkage measurement has to be done with the Length comparator with the . This version, conforming to UNI 6687, slightly differ from other version and it is used for the same determination.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003