


Product Name :
Hydraulic Test Pump

Product Code :
NLAB-MACHINELABWK28007



Description :

Hydraulic Test Pump

Technical Specification :

Hydraulic test pump designed for testing pressure cap. of different materials like drums, valves, boilers, barrels, tanks, for pressure vessels and cylinders complete with stand, pressure gauge, check valve and outlet valve with water tank having following specification.

Test Tank Capacity Kg. CM2

3000

2000

2000
3500
7000
10000

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003