

Product Name :
Hydraulic Valves and Fittings Test Stand

Product Code :
TK198



Description :

Hydraulic Valves and Fittings Test Stand

Technical Specification :

Test stand on which to mount industrial valves and fittings
Pressure testing of valves and fittings
Hand-operated piston pump to generate the test pressure, a return valve to relieve the system pressure, and a manometer for pressure measurement
2 different sizes of mounting flange with blank flange and flange seal
Connection of pump and test flange via pressure hose
Test medium: water
Mobile frame with collector tray and ball valve to drain
Water storage tank
Part of the Practice Line for assembly, maintenance and repair

Technical Data:

Piston pump with tank
test pressure: 0-60bar
tank capacity: 12L
manometer: 0-60bar
Mounting flanges for valves and fittings under test

Dimensions and weight:

Length x Width x Height: 1000x750x1200mm

Weight: 80kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003