

**Product Name :**  
Hydraulics Bench

**Product Code :**  
NLAB-VOCATIONALTVET19004



**Description :**

Hydraulics Bench

**Technical Specification :**

Water Circulation 1 HP Pump, Crompton.

Flow Measurement: Using Measuring Tank, Capacity 40 Ltrs.

Sump Tank Capacity

Stop Watch Electronic.

Control Panel Comprises of Standard make On-Off Switch, Mains Indicator, etc

Instruction Manual: An English instruction manual will be provided along with the Apparatus Tanks and Top Tray will be made of Stainless Steel

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

