

Product Name :
Hydraulics Bench

Product Code :
NLAB-VOCATIONALTVET19004



Description :

Hydraulics Bench

Technical Specification :

Water Circulation 1 HP Pump, Crompton.

Flow Measurement: Using Measuring Tank, Capacity 40 Ltrs.

Sump Tank Capacity

Stop Watch Electronic.

Control Panel Comprises of Standard make On-Off Switch, Mains Indicator, etc

Instruction Manual: An English instruction manual will be provided along with the Apparatus Tanks and Top Tray will be made of Stainless Steel

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

