

**Product Name :**  
Hydrostatic Test Bench

**Product Code :**  
NLAB-EDUCATIONTVET14003



**Description :**

Hydrostatic Test Bench

**Technical Specification :**

Bench intended to introduce the behavior of liquids under hydrostatic conditions (fluid at rest) and the application of these principles to pressure measurement using manometers.

Technical specifications:

Vertical tube with variable section, 460 mm scale

Demonstration of different types of pressure gauges:

Single tube piezometric pressure gauge

Inclined pressure gauge, angle of 5/30/60/90 °

Manometer with branches

U-tube manometer (air on liquid)

U-tube manometer (liquid on liquid)

Inverted U-tube manometer

Use of a vernier level meter

Dimensions (HxWxD): 1.25 x 0.425 x 0.46 m

Gross weight: 30 kg

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003