

Product Name :
Hydrostatic Weighing Frame

Product Code :
NLAB-EDUCATIONTVET11014



Description :

Hydrostatic Weighing Frame

Technical Specification :

Steel frame with crank for mounting the water tank

Supplied with plastic container 370x370x330 mm

Weighing basket

Hook for cubes and cylinders

Electronic platform scale

Capacity 15 kg

0.2g reading

Removable stainless steel platform 377 x 311 mm

4.3 "TFT WQVGA color graphics display

Application: Weighing, percentage weighing, piece counting, weight control,

Weighing hook under the scale

220V 50Hz power supply

Calibration mass 5kg, class F2, with certificate

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003