

Product Name :
Hydrostatic Weighing Frame

Product Code :
NLAB-EDUCATIONTVET13013



Description :

Hydrostatic Weighing Frame

Technical Specification :

Steel frame with crank for mounting the water tank

Supplied with plastic container 370x370x330 mm

Weighing basket

Electronic platform scale

Capacity 15 kg

0.2g reading

Removable stainless steel platform 377 x 311 mm

4.3 "TFT WQVGA color graphics display

Application: Weighing, percentage weighing, piece counting, weight control,

Weighing hook under the scale

220V 50Hz power supply

Calibration mass 5kg, class F2, with certificate

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003