

Product Name :
Ice Pack 0.3

Product Code :
NLABMEDICAL-M2310001



Description :

Ice Pack 0.3

Technical Specification :

Normal Capacity: 0.3 Litres

Volume of Ice Pack: 0.305 Litres

External Dimension (Approx.): 160 x 91 x 34 mm

Weight Empty: 76gm

Weight Filled: 363gm

Material: HDPE (High Density Polyethylene)

Performance:

Shipping Dimensions

Gross Weight (Approx.): 22kg

Size of carton (LxWxH) (approx.): 85 x 75 x 25cm

Volume of each carton (approx.): 0.16 CBM

No. of Units/20' Container: 41760

No. of Units/40' Container: 84960

No. of Units/40' High cube Container: 103680

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003