

Product Name :
ICU Bed Electric Remote Controlled

Product Code :
NLABMEDICAL-M2790002

Description :

ICU Bed Electric Remote Controlled India, ICU Bed Electric Remote Controlled manufacturer India, suppliers, Medical ICU Bed Electric Remote Controlled India, Hospital ICU Bed Electric Remote Controlled India

Technical Specification :

Base and mainframe work of precise steel tubes.

Base mounted on 125cm DIA

Swivelling castors, Central locking.

The backrest with CRP and knee rest as also trendelenberg

Reverse trendelenberg movement controlled by the Actuators

To achieve smooth and efficient operation without any jerks.

The level height of Bed can be adjusted to facilitate close

Observation by the Medical staff.

ABS bows easily removable for

Access to the patient.

Fitting: SS saline rod with four location,

protective railing

Finish in Epoxy powder coated.

Overall Size:

2160 L x 1015 W x 530 to 760 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003