

Product Name : ICU Bed Electric Remote Controlled	Product Code : NLABMEDICAL-M2790002
Description : ICU Bed Electric Remote Controlled India, ICU Bed Electric Remote Controlled manufacturer India, suppliers, Medical ICU Bed Electric Remote Controlled India, Hospital ICU Bed Electric Remote Controlled India	
Technical Specification : Base and mainframe work of precise steel tubes. Base mounted on 125cm DIA Swivelling castors, Central locking. The backrest with CRP and knee rest as also trendelenberg Reverse trendelenberg movement controlled by the Actuators To achieve smooth and efficient operation without any jerks. The level height of Bed can be adjusted to facilitate close Observation by the Medical staff. ABS bows easily removable for Access to the patient. Fitting: SS saline rod with four location, protective railing Finish in Epoxy powder coated. Overall Size: 2160 L x 1015 W x 530 to 760 mm	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003