

<b>Product Name :</b> Illuminated Zoom Microscope	<b>Product Code :</b> NLAB-EDUSCHOOLEQ173089
<b>Description :</b>	
Illuminated Zoom Microscope	
<b>Technical Specification :</b>	
<p>It is an ideal pocket microscope for Electronic Engineers, Hobbyists, Printing, Plate making, Scientists and Students.</p> <p>This is our smallest illuminated zoom microscope that provides 100 times linear magnification 10,000 times area magnification).</p> <p>Lenses are made of glass and provide clear images.</p> <p>Box measures 88mm x 40mm x 28mm. Each microscope comes with a synthetic leather pouch.</p>	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003