


Product Name : Illuminated Zoom Microscope	Product Code : NLAB-EDUSCHOOLEQ173089
--	---



Description :

Illuminated Zoom Microscope

Technical Specification :

It is an ideal pocket microscope for Electronic Engineers, Hobbyists, Printing, Plate making, Scientists and Students.

This is our smallest illuminated zoom microscope that provides 100 times linear magnification 10,000 times area magnification).

Lenses are made of glass and provide clear images.

Box measures 88mm x 40mm x 28mm. Each microscope comes with a synthetic leather pouch.

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003