

Product Name :
Impact Testing Machine

Product Code :
TN845



Description :

Impact Testing Machine

Technical Specification :

Impact Testing Machine It is designed for conducting Izod, Charpy test. The Pendulum is mounted on antifriction bearings. It has two starting positions, the upper one for Charpy and the lower one for Izod testing. Strike velocity of Pendulum:5.3m/sec for charpy 3.8m/sec for Izod, fixture for charpy and izod , Dial type.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003