

Product Name :
Incu Mixer MP Incubated Plate Vortexer

Product Code :
NLAB-EDUSCHOOLEQ45016



Description :

Incu Mixer MP Incubated Plate Vortexer

Technical Specification :

Incu Mixer MP is designed to provide an optimal method for vortexing and incubating microplates.

A compact, 2mm orbit ensures thorough mixing in microplates without splashing.

Simultaneous heating from both the bottom surface and the lid provide temperature uniformity throughout the sample.

Additionally, the heated lid minimizes evaporation and prevents condensation.

The Incu-Mixer MP series is ideal for a wide variety of molecular biology applications including immunochemical reactions, enzyme and protein analysis, and microarray analysis.

Designed for application flexibility, the Incu-Mixer series is available in two, and four-plate capacities, both capable of accepting all standard plates up to 40mm high.

Speed, time and temperature are all displayed in real time on the large LCD screen while an internal validation system constantly monitors all parameters for accuracy.

Independent time, temperature and speed controls allow the mixer to function as 3 machines in one; a multi-

plate incubator, vortexer, or a combination of the two.

<p>Ambier Range 70°C</p>	
<p>Orbit:</p>	
<p>Mixing Speed (4-Plate) 100 to 1200rpm (2-Plate)</p>	
<p>Temp. Increment:</p>	
<p>Temp. Accuracy:</p>	
<p>Temp. Uniformity:</p>	
<p>Speed Increment:</p>	
<p>Can Plate Height:</p>	

Dimensions (W x D x H) 11 x 10.6 x 5.5 in. / 28 x 27 x 14 cm (2-Plate)	
Weight (kg)	
Timer: to 99 hours 99min. / continuous	
Electrical 100 to 240V, 50-60Hz, 600W (4-Plate) 100 to 240V, 50-60Hz, 500W (2-Plate)	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003