

**Product Name :**  
Incu Shaker Mini with Flat Mat Platform

**Product Code :**  
NLAB-EDUSCHOOLEQ45013



**Description :**

Incu Shaker Mini with Flat Mat Platform

**Technical Specification :**

The Incu-Shaker Mini is the most compact shaking incubator in its class.

The standard platform (included) features a non-slip, rubber coated surface, ideal for tissue culture flasks, Petri dishes, and staining trays.

Despite the remarkably small footprint, a large 11.5 x 9.5 in. workspace accepts several platforms and a variety of flasks, tubes and other common vessels.

The unique, magnetic attachment method is the easiest way to instantly change between flask clamps of different sizes.

The convenient Magic Clamp platform is also available for use with Erlenmeyer flasks and test tube racks.

Spindle, 30 to 300 rpm	
------------------------	--

Orbit / 3/4 in.	
Temp. Range: 70°C	
Temp. Increment:	
Temp. Accuracy:	
Temp. Uniformity:	
Speed Increment:	
Platform Dimensions 29 cm	
Operating Temp. Range	
Dimensions (W x D x H) 41x 33 cm	

Weight: 11.5kg	
Height, 630mm Open:	
Timer: to 48 hours / continuous	
Electric: 230V, 50-60Hz, 300W	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003