

Product Name :
Indentation Hardness Tester Machines

Product Code :
NLAB-TESTINGLABM70001



Description :

Indentation Hardness Tester Machines

Technical Specification :

The feel of softness of any flexible foam is quantitatively defined by its indentation hardness. Indentation hardness of foam is determined by pressing a circular indenter at a specified speed against a block of foam and finding the force needed to compress it to a specified percentage of its initial thickness.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003