

Product Name : Induction Chines Wok	Product Code : NLAB-VOCATIONALTVET27099
---	---

--

Description :

Induction Chines Wok

Technical Specification :

Finishing Type : S.S.
Features : induction
Production Capacity (pcs per min.)
Power (W) : 3.5
Weight (kg) : 15 KG
Voltage (V) : 220V

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003