

<b>Product Name :</b> Induction Cooker	<b>Product Code :</b> NLAB-VOCATIONALTVET270101
<b>Description :</b> Induction Cooker	
<b>Technical Specification :</b>  Power : 3500 W  Material : Stainless Steel  Application : Pastry shops Delicatessen Supermarkets	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003