

Product Name :
Industrial Hydraulics Training Equipment

Product Code :
NLAB-VOCATIONALTVET19003



Description :

Industrial Hydraulics Training Equipment

Technical Specification :

Universal Testing Frame constructed of a twinned steel channel frame fastened with high tensile fixings.

This creates an overall structure size of 4.61m long, 2.53m high and 1.2m wide and a working space of 4.0m long x1.6m high.

The 0.60m spacing between the verticals allows for wide test specimens to be tested.

The base feet distribute the self-weight of about 1 tonne to four anti-vibration levelling feet.

The top channel members are used to carry either a single or double acting hydraulic ram Carriage.

All these accessories enable specimens to be set up quickly and easily in a variety of arrangements.

A set of hollow steel dial gauge support bars should also be included

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003