

**Product Name :**  
Industrial Refrigeration Trainer

**Product Code :**  
NLAB-TECHNICALAB11003



**Description :**

Industrial Refrigeration Trainer

**Technical Specification :**

Industrial Refrigeration Trainer consists of the following:

- 3 HP semi-hermetic compressor
- Refrigerant R-134A (other refrigerants available)
- Crankcase heater with automatic control system
- Oil circulating pump
- Two evaporators (forced-air type) with solenoids, timers, etc.
- Two water-cooled condensers (tube-in-tube and shell-in-tube)
- Gas bypass systems
- Crankcase pressure regulator with gauges

- Operating instructions, student experiment manual and teacher manual
- Text: Modern Refrigeration and Air Conditioning with study guide and lab manual

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003