

Product Name :
Infant Bed Child Cot S.S Adjustable Height

Product Code :
NLABMEDICAL-M3330005

Description :

Infant Bed Child Cot S.S Adjustable Height

Technical Specification :

Infant Bed / Cot.

Height adjustable approx. 850 - 1000 mm.

Completely made of stainless steel.

4 swivel castor (all with brakes).

Size: 750x400x850/1000mm

1000x400x850/1000mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003