

Product Name :
Infant Pediatric Cot

Product Code :
NLABMEDICAL-M2490010



Description :

Infant Pediatric Cot

Technical Specification :

Specification :-

Infant Pediatric Bed is manufactured with heavy round pipes.

Stainless Steel top with Mattress.

Body section is mounted on 10cm castors, 2 with Brakes.

Height adjustable approx. 85 - 100 cm

Product Size: 90 L x 50 W x 85-100 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Infant Pediatric Cot at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003