


<b>Product Name :</b> Infinity Corrected Digital Trinocular Microscope	<b>Product Code :</b> NLAB-EDUSCHOOLEQ1730110
---	--



**Description :**

Infinity Corrected Digital Trinocular Microscope

**Technical Specification :**

The digital microscope is equipped with a digital camera on a trinocular head, allowing users to go beyond conventional image processing functions to fully meet the requirements of your image analysis.

Equipped with a 3w LED lamp, field diaphragm, Abbe condensers and other optical components including a standard Kohler illumination system providing uniformly bright and free-from-glare specimen illumination which provides a high image quality with brighter and sharper resolutions for photomicrography.

The Infinity Corrected Microscope offers superior performance. Equipped with infinity corrected objectives, the built to withstand the demands of daily work in laboratory and classroom environments.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003