

Product Name :
Initial Surface Absorption Test (ISAT) Apparatus

Product Code :
NLAB -CONCRLAB120033



Description :

Initial Surface Absorption Test (ISAT) Apparatus

Technical Specification :

The rate of flow of water per unit area into a concrete surface when subjected to a constant head of 200mm is measured. The unit consists of a capillary tube mounted on a scale, a water reservoir & connecting tubes. Easy to use, mounted on a stand.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003