

Product Name :
Inside Micrometer

Product Code :
TN712



Description :

Inside Micrometer

Technical Specification :

Inside Micrometer
Carbide measuring faces
Equipped with a constant measuring-force device
An inside-length standard is required for accurately setting the micrometer
Range (5-30mm)
Resolution (0.001mm)
Accuracy 5 micron.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003