

<b>Product Name :</b> Instrument Cabinet	<b>Product Code :</b> NLABMEDICAL-M2730007
<b>Description :</b>	
Instrument Cabinet	
<b>Technical Specification :</b>	
Instrument cabinets are usually used in laboratories, hospital, research centers for keeping the drugs, medicines and lab chemicals safe	
These are manufactured in various designs such as with glass doors or without glass doors, with standard specifications or as per customer's specifications	
Frame Work Constructed with Heavy duty Ms Steel Sheet	
Two Glass Doors Available With Individual Key Locking	
Three Adjustable Glass Selves	
One Bottom Cabinet with Key Locking Facility	
Overall Size:- {1800(L)×1050(W)×450(D)} mm.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003