

**Product Name :**  
Insulating Oil Tester Fully Automatic

**Product Code :**  
NLAB-TESTINGLABM23003



**Description :**

Insulating Oil Tester Fully Automatic

**Technical Specification :**

Specification:

Fully automatic test cycle.

Automatic change of sample number of test.

Automatic calendar and time running even with power off.

Recall and print last 99 test results.

Programmable: Number of tests, settling time, stir time, waiting time.

LCD displays: Test status, HV rise, HV fall, BDV of previous test and test results.

Test voltage : 60, 80, 100 kV

Built in printer

Zero start safety and test chamber access interlock.

Ergonomic deesign, completely user friendly-just 36kg in weight.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003