

**Product Name :**  
Insulating Oil Testers

**Product Code :**  
NLAB-TESTINGLABM23005



**Description :**

Insulating Oil Testers

**Technical Specification :**

Self-contained

Compact and portable

Provide easy and accurate measurement of the breakdown voltage of insulating oil

Safe and easy operation

**Specifications:**

Three motorised rates of rise of voltage: 500/2000/3000 V per second

Stop/dwell rate of rise positions

Output voltages to 60 and 100 kV

Tests to: IEC 156, ASTM D1816, ASTM D877, JIS C2101, BS 5874, VDE 0370

Automatic HT shut down within 20 milliseconds

Digital memory kV meter holds breakdown voltage reading

Test voltage can be held at any value up to maximum

Zero start safety and test chamber access interlock

Fail Indicator Single unit construction with transparent window

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003