

**Product Name :**  
Interlocking Cubes

**Product Code :**  
NLAB-EDUSCHOOLEQ169024



**Description :**

Interlocking Cubes

**Technical Specification :**

1cm Interlocking Cubes

Perfect for counting, measuring and patterning. Plastic 1cm cubes in 10 colors come in a reusable Box.

Set of 1000 Pcs. in 10 color

2cm Interlocking Cubes

(Set of 64 Pcs. in 1 Color)

2cm Interlocking Cubes

(Set of 128 PCS. In 2 Color)

1cm Interlocking Cubes

Perfect for counting, measuring and patterning. Plastic 1 cm cubes in 10 colors come in a reusable Pouch. (Set of 1000 PCS, in 10 Color)

2cm Interlocking Cubes

Easy to connect on all sides, these interlocking 2 cm cubes come in 5 bright colors. Set of 625 Pcs. in 5 colors.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003