

Product Name :
Intermediate Vernier Microscope

Product Code :
NLAB-EDUSCHOOLEQ173038

Description :

Intermediate Vernier Microscope

Technical Specification :

The movement is controlled by a knob at the side of the base, and can be fixed in any position with a clamping screw, and a Vernier scale allows precise measurement to 0.01 mm.

The microscope itself features a 7x magnification Ramsden objective, 3x magnification 50 mm FL eyepiece, and a cross-piece graticule, with rack-and-pinion focusing.

The carriage has two locations for the microscope, which allows it to be mounted with its axis in any position relative to travel.

Supplied in a polished wooden case.

The unit is built on a heavy cast-iron base.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003