

Product Name :
Inverted Microscope

Product Code :
NLAB-EDUSCHOOLEQ1730102

Description :

Inverted Microscope

Technical Specification :

Precise fine focus controls move in 0.002 mm scale value increments. Phase Contrast models include a Centering Telescope.

The microscope comes with phase contrast objectives & condenser, 35mm and 65mm culture dish holders, and a trinocular photo port which allows the user to capture images with an add-on camera.

Featuring Binocular or Trinocular heads and Long Working Distance (LWD) objectives in your choice of Plan Achromat or Plan Phase Contrast.

The over-sized stage offers precise and smooth ball-bearing movement, for clear, steady viewing while scanning the field of view.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003