

Product Name :
IoT Trainer System

Product Code :
EC748



Description :

IoT Trainer System

Technical Specification :

IoT trainer system
IOT LAB Internet Thing Explore Trainer
Power 230V
Usage/Application: Laboratory
Type: Digital
Automation Grade: Automatic
Current: 2AMP
Color: MULTI
Frequency Range: 550HZS
IOT Base Smart Street Light Trainer
Design MULTI
Display Type: DIGITAL
IOT Sensor For Smart Environment Application Trainer
IOT Base Wireless Communication Module For Interfacing With Microcontroller Trainer

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments