

Product Name :
Isolation Gown

Product Code :
NLABMEDICAL-M3520001



Description :

Isolation Gown

Technical Specification :

Size: 115x127cm, 120x140cm Standard weight: 30-40 g/sq.m., Color: white, green, blue, yellow With or without poly-coating. Elastic cuffs or knitted cuffs. Attached waist ties and neck ties. Light and hygienic. Ideal for work in medical facilities.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003