

Product Name :
Isothermal Continuous Stirred Tank Reactor India

Product Code :
NLAB-VOCATIONALTVET29002



Description :

Isothermal Continuous Stirred Tank Reactor

Technical Specification :

Saponification Reaction study in Isothermal Continuous Stirred Tank Reactor (CSTR).

To determine the Reaction Rate Constant and effect of temperature on it.

Utilities Required:-

Electricity Supply: I Phase, 220VAC, 1.5kW.

Water Supply and Drain.

Instruments, Laboratory Glassware, and Chemicals required for analysis as per the system adopted.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com