

Product Name :
Isothermal Semi Batch Reactor India

Product Code :
NLAB-VOCATIONALTVET29004



Description :

Isothermal Semi Batch Reactor

Technical Specification :

Reactor: Material Stainless Steel, Capacity 3 Ltrs. (approx).

Water Bath: Material Stainless Steel, Double wall, insulated with Ceramic Wool.

Piping: Material Stainless Steel and PVC.

Pressure Regulator: 0-2 Kg/cm²

Pressure Gauge: Bourdon type 0-2 Kg/cm²

Stop Watch: Electronic

Heater: Nichrome wire Heater

Stirrer: Material Stainless Steel Impeller and shaft coupled with FHP motor

Feed Tank: Material Stainless Steel, Capacity - 20 Ltrs.

Feed Circulation: By compressed air.

Flow Measurement: Rotameter.

Control panel comprises of Digital Temp. Controller: 0-199.9°C

With Standard make on/off switch, Mains Indicator etc.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003