

**Product Name :**  
IV Stand Double Hook Metal Base 3 Wheels

**Product Code :**  
NLABMEDICAL-M2640003



**Description :**

IV Stand Double Hook Metal Base 3 Wheels

**Technical Specification :**

Specification :-

Frame Work of IV Stand is made of round tubes.

It comes with Base with 5 cms castors.

Base is made of Metal.

The adjustable IV rod is made of Stainless Steel  
with double hook.

The outer Pipe is made of Metal.

The height is adjustable from 137 to 240 Cms.

We are Manufacturer, Exporters and Suppliers of IV Stand Double Hook Metal Base 3 Wheels at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003