

Product Name :
IV Stand Double Hook Metal Base 5 Wheels

Product Code :
NLABMEDICAL-M2640001



Description :

IV Stand Double Hook Metal Base 5 Wheels

Technical Specification :

Frame Work of IV Stand is made of round tubes.

It comes with a Base with 5 cms castors.

The base is made of Metal.

The adjustable IV rod is made of Stainless Steel
with a double hook.

The outer Pipe is made of Metal.

The height is adjustable from 137 to 240 Cms.

We are Manufacturers, Exporters, and Suppliers of IV Stand Double Hook Metal Base 5 Wheels at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003