

**Product Name :**  
IV Stand Double Hook PU Base

**Product Code :**  
NLABMEDICAL-M2640002



**Description :**

IV Stand Double Hook PU Base

**Technical Specification :**

Specification :-

Frame Work of IV Stand is made of round tubes.

It comes with Base with 5 cms castors.

Base is made of PU (Polyurethane).

The adjustable IV rod is made of Stainless Steel with double hook.

The outer Pipe is made of Metal.

The height is adjustable from 137 to 240 Cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of IV Stand Double Hook PU Base at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003