

**Product Name :**  
IV Stand Four Foldable Hook Metal Base 3 Wheels

**Product Code :**  
NLABMEDICAL-M2640006



**Description :**

IV Stand Four Foldable Hook Metal Base 3 Wheels

**Technical Specification :**

Specification :-

Frame Work of IV Stand is made of round tubes.

It comes with Base with 5 cms castors.

Base is made of Metal.

The adjustable IV rod is made of Stainless Steel with foldable four hook.

The outer Pipe is made of Metal.

The height is adjustable from 137 to 240 Cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of IV Stand Four Foldable Hook Metal Base 3 Wheels at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003