

Product Name :
IV Stand Four Foldable Hook PU Base

Product Code :
NLABMEDICAL-M2640005



Description :

IV Stand Four Foldable Hook PU Base

Technical Specification :

Specification :-

Frame Work of IV Stand is made of round tubes.

It comes with Base having 5 cms castors.

Base is made of PU (Polyurethane).

The adjustable IV rod is made of Stainless Steel with foldable four hook.

The outer Pipe is made of Metal.

The height is adjustable from 137 to 240 Cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of IV Stand Four Foldable Hook PU Base at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003