

Product Name :
K-Slump Tester

Product Code :
NLAB -CONCRLAB120031



Description :

K-Slump Tester

Technical Specification :

The operation is very simple: insert the tester into the concrete up to the level of the disc, after 60 seconds, a measuring rod is lowered onto the surface of the concrete and the K-slump is read directly on a scale. The calibrated hollow tube has a dia. of 20 mm.

K-Slump Tester for in-place measurements indicates correlation to the slump test. The probe determines workability of concrete and the degree of compaction. Includes correlation chart and instructions.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003